

SMZ 10 Series

Ample Scientific




Features

- Zoom ratio of 1:10
- Crystal clear image and high three-dimensionality
- Rotatable head inclined at 45° for user's comfort
- Safe and durable guided-rail stand
- Trinocular head available

Specifications

Item	Specification
Head	Binocular (SMZ 10) or Trinocular (SMZ10 Trinocular)
Eyepiece	WF 10X, 20mm
Range of Magnification	1.9X—375X
Working Distance	31mm—320mm
Interpupillary Distance	50-75mm
Illumination	No illumination or Top & bottom LED illumination
Power	110V / 60Hz
Stand	Guided Rail Stand



Please contact Ample Scientific for available accessorizing options

Stereo Microscope

www.amplescientific.com