

SIM Series

Ample Scientific



Features

- Quality product at economical pricing
- Safe and durable stand with wall support
- Optional LED illumination on both top and bottom

Specifications

Item	Specification
Head	Binocular
Eyepiece	WF 10X, 19mm
Focus Distance	65 mm
Working Distance	95 mm
Interpupillary Distance	55-75 mm
Size	62 x 39.5 x 48 cm
Weight	4.2 Kg
Illumination	No illumination or Top & bottom LED illumination
Power	110V / 60Hz
Objectives	1X / 3X or 2X / 4X



Please contact Ample Scientific for available accessorizing options

Objective Magnification	Visual Field	Total Magnification
1X	19.0 mm	10X
2X	9.5 mm	20X
3X	6.3 mm	30X
4X	4.8 mm	40X

Stereo Microscope

www.amplescientific.com